

# Playground

Exported from [IcePanel](#) at 2024-02-22 21:55:00

## Table of contents


- [1. Context diagrams](#)
  - [1a. Monitor Me - Arch Kata](#)
- [2. Areas](#)
  - [2.1. MonitorMed](#)
- [3. Actors](#)
  - [3.1. Medical Professional](#)
  - [3.2. Nurse](#)
  - [3.3. Patient](#)
- [4. Systems](#)
  - [4.1. Monitoring Device](#)
  - [4.2. MonitorMe](#)
    - [4.2.1. App diagrams](#)
      - [4.2.1a. MonitorMe - App Diagram](#)
      - [4.2.1b. Nurse assigns medical professional for alerting purposes](#)
      - [4.2.1c. Nurse monitors vitals](#)
      - [4.2.1d. Medical Professional and/or nurse receive alert\(s\)](#)
      - [4.2.1e. MonitorMe device\(s\) monitor patient](#)
      - [4.2.1f. Medical Professional and/or nurse pulls snapshot of patient vitals](#)
      - [4.2.1g. Nurse onboards patient with devices](#)
      - [4.2.1h. Medical Professional and/or nurse configures alerts](#)
    - [4.2.2. Apps](#)
      - [4.2.2.1. Alerts](#)
      - [4.2.2.2. API Gateway](#)
      - [4.2.2.3. Devices](#)

- [4.2.2.4. MonitorMe Event Broker](#)
    - [4.2.2.5. Patients](#)
      - [4.2.2.5.1. Component diagrams](#)
    - [4.2.2.6. Raw Vitals](#)
    - [4.2.2.7. Snapshots](#)
    - [4.2.2.8. Subscriptions](#)
    - [4.2.2.9. Vital Events](#)
  - [4.2.3. Stores](#)
    - [4.2.3.1. MonitorMe \(Database\)](#)
      - [4.2.3.1.1. Component diagrams](#)
        - [4.2.3.1.1a. Diagram 1](#)
      - [4.2.3.1.2. Components](#)
        - [4.2.3.1.2a. Alert Rules \(Table\)](#)
        - [4.2.3.1.2b. Alerts \(Table\)](#)
        - [4.2.3.1.2c. Devices \(Table\)](#)
        - [4.2.3.1.2d. Patients \(Table\)](#)
        - [4.2.3.1.2e. Patients Alert Rules \(Table\)](#)
        - [4.2.3.1.2f. Patients Device \(Table\)](#)
        - [4.2.3.1.2g. Snapshots](#)
        - [4.2.3.1.2h. Subscriptions \(Table\)](#)
        - [4.2.3.1.2i. Users \(Table\)](#)
    - [4.2.3.2. MonitorMe Raw Vitals \(Database\)](#)
      - [4.2.3.2.1. Component diagrams](#)
        - [4.2.3.2.1a. MonitorMe Raw Vitals \(Database\) Component Diagram](#)
      - [4.2.3.2.2. Components](#)
        - [4.2.3.2.2a. Raw Vitals \(Table\)](#)
- [4.3. MonitorMe Device](#)
  - [4.3.1. App diagrams](#)

- [4.4. MyMedicalData](#)
- [4.5. Push Platform](#)
- [4.6. StayHealthy Mobile Devices](#)
  - [4.6.1. App diagrams](#)
    - [4.6.1a. Diagram 1](#)
  - [4.6.2. Apps](#)
    - [4.6.2.1. StayHealthy Android](#)
    - [4.6.2.2. StayHealthy iOS](#)
- [4.7. StayHealthy Monitoring Station](#)
  - [4.7.1. App diagrams](#)
    - [4.7.1a. StayHealthy Monitoring Station App Diagram](#)
  - [4.7.2. Apps](#)
    - [4.7.2.1. Monitoring Station App](#)

# 1. Context diagrams

## 1a. Monitor Me - Arch Kata



## 2. Areas

### 2.1. MonitorMed

No description

## 3. Actors

### 3.1. Medical Professional

No description

### 3.2. Nurse

No description

### 3.3. Patient

No description

# 4. Systems

## 4.1. Monitoring Device


No description

## 4.2. MonitorMe


No description

### 4.2.1. App diagrams

#### 4.2.1a. MonitorMe - App Diagram



#### 4.2.1b. Nurse assigns medical professional for alerting purposes




- Step 1: Nurse -> Monitoring Station App: Nurse selects medical professional in system to be alerted via Interacts with monitoring station app
- Step 2: Monitoring Station App -> API Gateway: App sends request to create a patient notification request via Create, manage patients, devices, alerts
- Step 3: API Gateway -> Subscriptions: Gateway sends request to subscriptions service via Create, accept request for patient alerts
- Step 4: Subscriptions -> MonitorMe (Database): Subscriptions service create subscription request for medical professional via CRUD subscriptions
- Step 5: Medical Professional -> StayHealthy iOS: Medical professional logs into the app via Views data and alerts
- Step 6: Medical Professional -> StayHealthy iOS: Medical professional accepts patient to monitor via Views data and alerts
- Step 7: StayHealthy iOS -> API Gateway: App sends push token along with the subscription request to system via Create, manage patients, devices, alerts
- Step 8: API Gateway -> Subscriptions: Gateway sends request to subscriptions via Create, accept request for patient alerts




- Step 9: Subscriptions -> MonitorMe (Database): Subscriptions microservice saves push token in subscriptions table via CRUD subscriptions

## 4.2.1c. Nurse monitors vitals



- Step 1: Nurse -> Monitoring Station App: Nurse monitors vitals via Interacts with monitoring station app
- Step 2: Monitoring Station App -> API Gateway: Sends request to see vitals via Create, manage patients, devices, alerts
- Step 3: API Gateway -> Vital Events: Vitals creates server sent events for each vital via Vital server sent events
- Step 4: API Gateway -> Monitoring Station App: App receives events and refreshes display via Create, manage patients, devices, alerts
- Step 5: Monitoring Station App: App determines patient's devices
- Step 6: Monitoring Station App -> Nurse: App shows nurse vitals for patient via Interacts with monitoring station app


#### 4.2.1d. Medical Professional and/or nurse receive alert(s)



- Step 1: MonitorMe Device -> MonitorMe Event Broker: Device measures value and pushes to system via Push data
- Step 2: MonitorMe Event Broker -> Alerts: Broker sends vital to alerts via Notify alerts of all vitals to be evaluated
- Step 3: Alerts: Alert checks rules
- Step 4: Alerts: Alert evaluates vital against rules
- Step 5: Alerts -> MonitorMe (Database): Alert saves alert in database via Read subscriptions (push tokens), read patient, read patient\_device, read device, crud alert rules, crud patient\_alert\_rules
- Step 6: Alerts -> MonitorMe (Database): Alert checks subscriptions via Read subscriptions (push tokens), read patient, read patient\_device, read device, crud alert rules, crud patient\_alert\_rules
- Step 7: Alerts -> Push Platform: Alert sends alert via Send alerts to be pushed
- Step 8: Push Platform -> StayHealthy Mobile Devices: Push platform sends notification to device via Push notifications


- Step 9: StayHealthy iOS: iOS alerts medical professional

#### 4.2.1e. MonitorMe device(s) monitor patient




- Step 1: MonitorMe Device -> MonitorMe Event Broker: Sends vitals to broker via Push data
- Parallel path: Send vitals to consumers
  - Step 1: MonitorMe Event Broker -> Raw Vitals: Broker sends vitals to Raw Vital service via Persist vitals
  - Step 2: MonitorMe Event Broker -> Alerts: Broker sends vitals to Alerts service via Notify alerts of all vitals to be evaluated
  - Step 3: MonitorMe Event Broker -> Vital Events: Broker streams vitals to Vital Events via Stream events to be sent back to frontends

#### 4.2.1f. Medical Professional and/or nurse pulls snapshot of patient vitals




- Step 1: Nurse: Nurse collects patient information.
- Step 2: Nurse -> Monitoring Station App: Initiates snapshot creation. via Interacts with monitoring station app
- Step 3: Monitoring Station App -> API Gateway: Makes the snapshot creation request. via Create, manage patients, devices, alerts
- Step 4: API Gateway -> Snapshots: Request is sent to Snapshots service. via Create snapshot
- Step 5: Snapshots -> MonitorMe Raw Vitals (Database): Patient's raw vitals are queried. via Query time series data based on device id and time
- Step 6: Snapshots -> MonitorMe (Database): Snapshot event is saved to the database. via Write, Read snapshots, Read patients, Read devices,
- Step 7: Snapshots -> MyMedicalData: Snapshot is sent to MyMedicalData. via Send snapshots

#### 4.2.1g. Nurse onboards patient with devices



- Step 1: Nurse: Nurse collects patient information
- Step 2: Nurse: Nurse prepares devices for patient
- Step 3: Nurse -> Monitoring Station App: Nurse registers patient with devices on the Monitoring Station App via Interacts with monitoring station app
- Step 4: Monitoring Station App -> API Gateway: App sends patient and device registration details to MonitorMe via Create, manage patients, devices, alerts
- Step 5: API Gateway -> Patients: Gateway sends request to Patients microservice via Create, manage patients
- Step 6: Patients -> MonitorMe (Database): Patients microservice creates patient, creates new device if not existing, relates patient to device via Write, Read patients

#### 4.2.1h. Medical Professional and/or nurse configures alerts



- Step 1: Nurse -> Monitoring Station App: Nurse selects patient to configure alerts for via Interacts with monitoring station app
- Step 2: Nurse -> Monitoring Station App: Nurse configures alerts based on vital signs via Interacts with monitoring station app
- Step 3: Monitoring Station App -> API Gateway: App sends request via Create, manage patients, devices, alerts
- Step 4: API Gateway -> Alerts: Gateway sends request to alerts microservice via Create, manage alert rules
- Step 5: Alerts -> MonitorMe (Database): Alert service save new alert rules via Read subscriptions (push tokens), read patient, read patient\_device, read device, crud alert rules, crud patient\_alert\_rules
- Step 6: Alerts -> MonitorMe (Database): Alert service associates alert rule with patient via Read subscriptions (push tokens), read patient, read patient\_device, read device, crud alert rules, crud patient\_alert\_rules
- Step 7: Alerts: Alert rules are synced in memory

## 4.2.2. Apps

### 4.2.2.1. Alerts

Determines alerts based off of alert rules and sends alerts to subscribers

Golang

- Create alert rules given patient
- Send alerts to push platform
- 

### 4.2.2.2. API Gateway

NGINX

No description

### 4.2.2.3. Devices

MonitorMe devices (responsible for registering devices, associating devices, and deleting devices)

Golang

- Registers devices
- Associate device to patient
- Delete device

### 4.2.2.4. MonitorMe Event Broker

Apache Kafka

No description

### 4.2.2.5. Patients

Create, manage patients

Golang

- Create new patient

#### 4.2.2.5.1. Component diagrams

### 4.2.2.6. Raw Vitals

Golang

No description

### 4.2.2.7. Snapshots

Golang

- Create snapshot for patient
  - Look up patient and their devices
  - Look up the patient's usage of device timeframe
  - Look up raw vitals
  - Package and send data

## 4.2.2.8. Subscriptions

Manage, accept patient notification (alerts) requests

Golang

- Create request for patient notifications
- Accept medical professional request for patient notifications
  - Receives push token
- 

## 4.2.2.9. Vital Events

Golang

- Reads raw vitals back to applications
  - Can be filtered based on vitals and patient/device ID
- 

## 4.2.3. Stores

### 4.2.3.1. MonitorMe (Database)


PostgreSQL

No description

#### 4.2.3.1.1. Component diagrams

##### 4.2.3.1.1a. Diagram 1





## 4.2.3.1.2. Components

### 4.2.3.1.2a. Alert Rules (Table)

No description

### 4.2.3.1.2b. Alerts (Table)

No description

### 4.2.3.1.2c. Devices (Table)

No description

### 4.2.3.1.2d. Patients (Table)

No description

### 4.2.3.1.2e. Patients\_Alert\_Rules (Table)

No description

### 4.2.3.1.2f. Patients\_Device (Table)

No description

#### 4.2.3.1.2g. Snapshots

No description

#### 4.2.3.1.2h. Subscriptions (Table)

No description

#### 4.2.3.1.2i. Users (Table)

No description


### 4.2.3.2. MonitorMe Raw Vitals (Database)

PostgreSQL

No description

#### 4.2.3.2.1. Component diagrams

##### 4.2.3.2.1a. MonitorMe Raw Vitals (Database) Component Diagram



#### 4.2.3.2.2. Components

#### 4.2.3.2.2a. Raw Vitals (Table)

- Device ID
- Timestamp
- Device Type (Vital type)
- Value

### 4.3. MonitorMe Device

External

8-Patient Monitoring Equipment

#### 4.3.1. App diagrams

### 4.4. MyMedicalData

External

No description

### 4.5. Push Platform

External


Apple Push Notifications and/or Firebase

### 4.6. StayHealthy Mobile Devices

No description

#### 4.6.1. App diagrams

##### 4.6.1a. Diagram 1



## 4.6.2. Apps

### 4.6.2.1. StayHealthy Android

React Native

No description

### 4.6.2.2. StayHealthy iOS

React Native

No description

## 4.7. StayHealthy Monitoring Station

Monitoring screen in nurses station

### 4.7.1. App diagrams

#### 4.7.1a. StayHealthy Monitoring Station App Diagram

## StayHealthy Monitoring Station



### Monitoring Station App

[App: React Native]

Monitoring screen in nurses station

[System]

## 4.7.2. Apps

### 4.7.2.1. Monitoring Station App

React Native

No description